


# Product Information Guide

**widex mind**  
the pleasure of hearing

- Dual Integrated Signal Processing Technology with Two Layers: Listening and Awareness
- 15 Band TruSound Processing
- 15 Band High Definition Locator Microphone
- Multidirectional Feedback Canceling
- Zen Relaxation Program
- Speech Enhancer
- Up to 5 Customized Listening Programs
- Smart Speak Signaling Mechanism
- Audibility Extender
- In-Situ RECD (*m4-9 model only*)
- Enhanced Sound Diary II
- ClearBand Receiver (*m4-m-CB model only*)
- Dedicated RC4 Remote Control Option
- LED Indicator (*m4-19 model only*)

## m4-CIC

15 Channel, Omnidirectional




Styles Available	CIC, Mini Canal
Peak Output	112 dB SPL
Max Gain*	51 dB
Battery Size	10
Battery Life***	125 hours
Battery Drain	.75 mA
Telecoil TLS	N/A
Volume Control	Remote Control Only
Mobile Phone Immunity (Urating)****	M3

	500	1000	2000	4000 Hz
Max HTL	65	75	75	75 dB

## m4-X

15 Channel, High Definition Locator




Styles Available	Canal, Half Concha, Full Shell
Peak Output	114 dB SPL
Max Gain*	58 dB
Battery Size	312
Battery Life***	160 hours
Battery Drain	.9 mA
Telecoil Output (Average) **	86 dB SPL
Volume Control	Optional
Mobile Phone Immunity (Urating)****	M3/T3

	500	1000	2000	4000 Hz
Max HTL	85	85	85	85 dB

## m4-9

Mini BTE, 15 Channel, High Definition Locator




Styles Available	Mini BTE
Peak Output	124 dB SPL
Max Gain*	57 dB
Battery Size	13
Battery Life***	315 hours
Battery Drain	.85 mA
Telecoil Output (Average) **	89 dB SPL
Volume Control	Optional
Mobile Phone Immunity (Urating)****	M3/T3

	500	1000	2000	4000 Hz
Max HTL	80	80	80	80 dB

## m4-19

Power BTE, 15 Channel, High Definition Locator




Styles Available	BTE
Peak Output	131 dB SPL
Max Gain*	70 dB
Battery Size	13
Battery Life***	325 hours
Battery Drain	.85 mA
Telecoil Output (Average) **	97 dB SPL
Volume Control	Optional
Mobile Phone Immunity (Urating)****	M3/T3

	500	1000	2000	4000 Hz
Max HTL	90	95	100	100 dB

## m4-m

Micro Size BTE, 15 Channel, High Definition Locator




Styles Available	Micro BTE
Peak Output	126 dB SPL
Max Gain*	58 dB
Battery Size	10
Battery Life***	110 hours
Battery Drain	.8 mA
Telecoil Output (Average) **	N/A
Volume Control	Remote Control Only
Mobile Phone Immunity (Urating)****	M3

	500	1000	2000	4000 Hz
Max HTL	75	80	80	80 dB

## m4-m CB

Micro size BTE, 15 Channel, High Definition Locator, ClearBand Receiver



Styles Available	Micro BTE
Peak Output	122 dB SPL
Max Gain*	55 dB
Battery Size	10
Battery Life***	105 hours
Battery Drain	.85 mA
Telecoil Output (Average) **	N/A
Volume Control	Remote Control Only
Mobile Phone Immunity (Urating)****	M3

	500	1000	2000	4000 Hz
Max HTL	75	80	80	80 dB

**Notes:**  
\*Real max gain depends on vent size and other variables such as the feedback canceling mode selected. \*\*Input 10 mA/m. \*\*\*Lifetime of the battery in test mode with ANSI 65dB/1kHz signal. \*\*\*\*Microphone (M) and telephone (T) immunity test ratings reported are for both 800-950 MHz and 1600 MHz – 2.4 GHz bands unless otherwise indicated. Immunity rating depends on the frequency band used by the wireless service provider. Contact your wireless service provider regarding this and any additional information or necessary phone settings for optimum compatibility.

